

ADVALVE

Bioprosthesis

Let Life Endure.

Technical Brochure



Advanced Cardiac Surgery

advancing the science of cardiothoracic surgery




Beyond Endurance, Beyond Excellence.

Spatial Modeling Technology

Composite Metal Polymer Stent ^{1,2}


Redistributes significant load to the polymer frame.
The Nitinol wire circuit acts as a vibration dampener.



Advanced Leaflet Technology

Advanced Digital Fabrication

Pericardial thickness detection ensures uniform prosthesis quality.
Even stress distribution across leaflet surface.



Tissue Fixation Technology

eGDGE Cross-Linking Minimizes Calcification.






Fig.1 - Accumulation of calcium in biomaterial implanted in animals (accelerated calcification model). Next to the name of the fixative is a diagram of its molecule, where the hydrocarbon chain; ● reaction groups —C—C—. At point "0" is the amount of calcium in the non-implanted biomaterial.³

Lower Calcification Risk: eGDGE vs GA

Experimental models show a pronounced inhibition of calcification with eGDGE.⁴



Biomechanical Characteristics³:

Enhanced Elasticity and Durability, Resistant to Long-Term Degradation.





Fig.2 - PV - Porcine Aortic Valve Leaflets, P - Pericardium, C- Cattle Arteries³

Anti-Bacterial Treatment:

Fibrous film of chlorhexidine bigluconate on the surface of an eGDGE-treated prosthesis.



Transprothesis Gradient





Fig.3 - Transprothesis gradient of ADVALVE prostheses treated with diepoxy compound (EGDGE) or glutaraldehyde (GA) during bench tests in the mitral position. *- significant differences (p<0.05) between the groups of GA and EGDGE bioprostheses.³



Minute Volume




Fig.5 - Minute volume of ADVALVE prostheses treated with diepoxy compound (EGDGE) or glutaraldehyde (GA), mitral position, bench tests.³

Clinical Outcomes²

Absence of Primary Tissue Failure (PTF)
EGDGE-treated bioprostheses vs GA-treated counterparts.




Fig.4 -Actuarial results of primary tissue failure (PTF) absence of bioprostheses in patients younger than 50 years old (according to various studies). KCC is a group of patients with epoxy-treated bioprostheses (own data) compared to Santini F. et al., Yu H. et al. Jamieson W. et al.²

Lower Immunogenicity^{2,3}

Compared to GA, EGDGE-treated tissues lower immune response and deterioration risk.

Compatibility with Advanced Modifications^{2,3}

The absence of excessive free aldehyde groups in EGDGE-treated biomaterials allow antithrombotic and antibacterial coatings without compromising integrity.

Advanced Anti-Calcification²



Effective Inhibition of Calcification

Diphosphonates inhibit calcium phosphate and hydroxyapatite crystallization.

Long-Term Protection

Prolonged Anti-calcification due to covalent binding of diphosphonates with EGDGE-treated groups, beneficial for paediatric & young patients.

Antibacterial Protection

Chlorhexidine immobilization resists gram-negative and gram-positive microbes & Improves valve replacement outcomes in Infectious endocarditis.


References

1. Karaskov, A. M., et al. "Uniline Next-Generation Russian Biological Prosthesis for Mitral Valve Replacement: First Experience." Paper presented at: Russian Journal of Thoracic and Cardiovascular Surgery Conference, 1-6 Feb. 2017, Russia.
2. Barbarash, L. S., and I. Y. Zhuravleva. "Evolution of Heart Valve Bioprostheses: Achievements and Challenges of Two Decades." Paper presented at: Scientific Research Institute of Complex Problems of Cardiovascular Diseases, 2012, Kemerovo, Russia.
3. Data on file.

** ADVALVE is manufactured by AMS under a technology license from Neocor for Uniline.

Technical Specifications

AORTIC ADVALVE Aortic



Nominal Specifications (mm)

AVA19s

AVA21s

AVA23s

AVA25s

D - Stent Diameter	19	21	23	25
D1 - External Sewing Ring Diameter	24.5	27	28	30
H - Total Height	16	16.5	17.0	17.5
H1 - Aortic Protrusion	13	13.5	14.5	15.0

MITRAL ADVALVE Mitral



Nominal Specifications (mm)

AVM26s

AVM28s

AVM30s

AVM32s

D - Stent Diameter	26	28	30	32
D1 - External Sewing Ring Diameter	38.0	40.0	42.0	43.5
H - Total Height	15.0	16.0	16.5	17.5
H1 - Ventricular Protrusion	12	13.5	15	15

